




# Why your choice of breast shield is key to successful pumping

## Understanding your breast anatomy

Breasts come in all shapes and sizes, and can change throughout your breastfeeding journey.<sup>1</sup> The size of your breast might not correspond to your nipple size, and each of your nipples can be different sizes too.




## What happens when you breastfeed or pump



Your milk ducts increase in size by 68% during a let down.<sup>2</sup> This is to accommodate all the milk flowing through them towards your nipple.



Your nipple diameter may increase temporarily by 2 to 3 mm.<sup>3</sup> So your breast shield tunnel needs to be slightly bigger than your nipple.



Your milk ducts lie close to the skin's surface. Pressing the breast shield too hard or using a tunnel that's too small can obstruct milk flow.<sup>4</sup>

## Why comfort makes a difference to successful pumping


**Stress and discomfort can hinder the hormone oxytocin, essential for the release of breast milk.<sup>5</sup> See our tips for more efficient pumping:**



Using a correctly fitting breast shield is important for you to be relaxed and comfortable while pumping and helps your milk to flow.<sup>6</sup>



Pumping shouldn't hurt. It is important to adjust the vacuum to the highest comfortable level during the expression phase to help remove more milk in less time.<sup>7</sup>



Your nipple should be centred and moving freely in the tunnel during pumping. If it doesn't, you need to choose another size breast shield.

**4.** Selecting a pump with overflow protection – sometimes called a closed system – will allow you to pump in a position that works for you.



Find our breast shield fitting guide at [www.medela.com/fittingguide](http://www.medela.com/fittingguide) and read more at [www.medela.com/pumping](http://www.medela.com/pumping)

**References:** 1. Cox DB et al. Exp Physiol. 1999;84:421-434. 2. Ramsay DT et al. Pediatrics. 2004;113:361-367. 3. Geddes DT et al. Early Hum Develop. 2008 Jul 1;84:471-477. 4. Geddes, DT. J Midwifery Womens Health. 2007;52:556-563. 5. Newton M, Newton NR. J Pediatr. 1948;33:698-704. 6. Jones E, Hillon S. J Neonatal Nurs. 2009;15;15:14-17. 7. Kent JC et al. Breastfeed Med. 2008;3:11-19.